INDIGO[®]

INSTRUCTION MANUAL

NH-81-00003 Rev C

CONTENTS

Introduction	1	Fire Safety	20
Terms and Definitions	1	Burn Prevention	20
Overview	3	Child Safety	21
Components	4	Cleaning Safety	21
Vaccine Carrier Use	5	Components	22
Components	5	Overview	22
Bag	6	Display	22
Cooler Lid	7	Charging and USB Ports	22
Display	8	Installation	23
Preparation for Use	9	Shipped Components and Tools	23
Normal Operation	9	Installation Instructions	23
Intermittent Use	12	Power Requirements	28
Bag Strap Adjustment	13	Operation	29
Data Download	14	Data Download	36
Battery Replacement	15	Cleaning	37
Cleaning	16	Signs of Damage	37
Signs of Damage	16	Troubleshooting	38
Troubleshooting	17		
Charging	19	Disposal and Recycling Guidelines	39
Safety Precautions	19		
Proper Installation and Maintenance	20	Warranty Terms and Conditions	40

Terms and Definitions

- **Absolute Time**

Coordinated Universal Time (UTC) time derived from an independent verified source, such as a cellular tower, GPS, or Internet time server.

- **Ambient Temperature**

Air temperature of the environment where equipment is stored.

- **Cold Life**

Estimated remaining time that the vaccine storage compartment will stay within 2 to 8 °C before the device will need to be recharged.

- **Cool-down**

The time required for the device to cool, starting when the vaccine carrier is switched on and ending when temperatures at all areas in the vaccine storage compartment are between 2 and 8 °C.

- **Equipment Monitoring Device**

A device that retrieves data from an appliance's logger and stores, analyzes, and communicates data.

- **Immunization**

The process whereby a person is made immune or resistant to an infectious disease, typically by the administration of a vaccine.

- **Outreach**

The practice of health workers delivering and administering vaccines in communities, typically those in remote locations, to improve vaccine access.

- **Relative Time**

A timestamp using a constant reference point, such as the time of device commissioning.

- **Secondary Packaging**

Packaging, such as a carton, containing several individual vaccine containers.

- **USB Mass Storage Device**

A USB device that follows protocols defined by the USB Implementors Forum to enable file transfer between the host and USB device.

- **Vaccine Storage Compartment**

The chamber in the vaccine carrier that is suitable for storing vaccines.

INDIGO SYSTEM

Overview

The Indigo system is designed to enable ice-free, freeze-free immunization outreach. Indigo utilizes ColdLoop™ technology to ensure five days of freeze-free cooling capacity in even the most extreme environments. This means the carrier stays at a vaccine safe temperature—between 2 and 8 °C—throughout five days of continuous operation at 43 °C ambient temperature. The portable vaccine carrier presents opportunities to extend capabilities of vaccinators and to immunize children in the most hard-to-reach communities. The system is designed primarily for use in conducting immunization outreach activities, although it may be used for other transportable refrigeration applications.

Components

Indigo Vaccine Carrier: Cooler, cooler lid, and bag are used together during immunization outreach.

- **Cooler:** The container used to refrigerate vaccines. The cooler integrates evaporative cooling technology with vacuum insulated stainless steel to maintain a vaccine storage compartment temperature of 2 to 8 °C for five days of continuous use at 43 °C ambient temperature. A cooler's fully or partially depleted cold life can be recharged using the Indigo Charger.
- **Cooler Lid:** The insulating lid on top of the cooler, which enables vaccine access,

The Indigo system has a unique set of features and capabilities that address the limitations of traditional ice-based carriers:


- **Extended cooling:** Provides vaccine-safe storage temperatures for five days of continuous use at 43 °C ambient temperature, without external power.
- **On/off switch:** The cold life of a charged vaccine carrier may be stored until needed.
- **Freeze-free:** Prevents vaccine exposure to freezing temperatures without action from the user.
- **Reusable and durable:** The vaccine carrier is rechargeable and made from durable materials.

displays temperature, shows estimated remaining cold life, and logs data.

- **Bag:** The case used for carrying and protecting the cooler, cooler lid, and other immunization supplies. The bag can be carried comfortably as a backpack or mounted on a vehicle.


Indigo Charger: Used to recharge the cooler after the cold life has been fully or partially depleted. The charger restores the cooling power of the cooler using a heating and cooling process and resets the cold life indicator on the cooler lid. A typical charge takes four to five hours.

Indigo System

Indigo Charger

Indigo Vaccine Carrier

Bag

Cooler Lid



Cooler

VACCINE CARRIER USE

Components

- 1) Cooler Inner Lid
- 2) Cooler Outer Lid
- 3) Display
- 4) Bag
- 5) Vaccine Storage Compartment
- 6) Cold Key
- 7) Valve


Bag

- 1) Upper Zipper:** Opened to access vaccines and view the cooler lid display.
- 2) Lower Zipper:** Used to remove the cooler and cooler lid for charging and maintenance.
- 3) Top Pocket:** Stores the foam vial holder.
- 4) Front Pocket:** Stores the USB On-The-Go

(OTG) adaptor for data download (see "Data Download"). May be used to store other items, such as immunization records and/or compatible data devices.

- 5) Strap Adjusters:** Used to adjust the strap length (see "Bag Strap Adjustment").


- 6) Tie-down Straps:** Used to strap securely to bicycle, motorbike, or vehicle.

Cooler Lid

NOTICE: The cooler lid has a battery intended to provide power for approximately two years and does not need to be recharged. See “Battery Replacement” section for instruction on replacing the battery.

- 1) **Cooler Inner Lid:** Removable to access vaccines.
- 2) **Cooler Outer Lid:** Remains on the cooler during immunization activities and is removable for charging and maintenance.
- 3) **Electronics Enclosure:** Houses lid electronics and extra battery (see “Battery Replacement”) and displays vaccine storage compartment temperature and other information.
- 4) **Temperature Probe:** Measures temperature of vaccine storage compartment.
- 5) **Cold Key:** Used to turn cooling on and off.
- 6) **Tether:** Secures the cooler outer lid to the cooler.
- 7) **USB-A Cable:** Connects to the charger during charging and connects to compatible data devices for data download during or after use (see “Data Download”).

Display


The cooler lid display will turn on when the cooler is turned on. The display will turn off when the cooler is turned off but can be

activated by opening or closing the cooler inner lid or plugging the USB-A cable into a compatible device.

- 1) **Temperature:** Allows a user to monitor the vaccine storage compartment temperature in degrees Celsius.
- 2) **Cold Life Indicator:** Displays the estimated remaining cold life. The number inside the circle shows the day of use, counting down from "5 DAYS" to "0 (zero) DAYS". The outer circle shows the portion remaining of the current day of use, with each segment representing 2.4 hours. After a full charge, estimated cold life will update to "5 DAYS" and the outer circle will show all segments.
- 3) **ON Indicator:** Appears if cold key is on the cooler valve.
- 4) **HIGH / LOW Alarms:** Displays alarm status when the cold key is on the cooler valve and the vaccine storage

compartment temperature exceeds 10 °C for 10 hours or drops below -0.5 °C for 60 minutes. To reset alarms, briefly open and close the inner lid, when vaccine chamber temperature is between 2 and 8 °C.


- 5) **USB Connected:** Appears when the cooler lid USB is successfully connected to the charger or a compatible data device (see "Data Download").
- 6) **OPEN Indicator / Alarm:** OPEN indicator, "OPEN", appears when the cooler inner lid is open. OPEN alarm symbol, "!" appears when the cooler inner lid is open for more than 60 minutes. Turns off when cooler inner lid is closed again.
- 7) **Replace Battery:** Appears when battery in the cooler lid reaches a low voltage and should be replaced (see "Battery Replacement").

Preparation for Use

Charge the cooler prior to initial use. A fully charged cooler will display “5 DAYS” of cold life on the cooler lid. See charger label and “Charging” manual section for instructions. Charge using the Indigo Charger only.

- 1) **Charger Label**
- 2) **Indigo Manual**

Normal Operation


NOTICE: Do not switch the cooler on prior to the intended use. Additional cool-downs from ambient temperature will deplete the cooler’s cold life.

NOTICE: This device is not intended for extended vaccine storage longer than five days at 43 °C ambient temperature.

- 1) **Check cold life indicator for cooler charge status:** Open the cooler inner lid to display the cold life indicator and check the cooler’s charge. A fully charged cooler will show “5 DAYS” of cold life with a complete circle around the number. A partially charged cooler will show a number less than five and/or a partial circle around the number.


- 2) Turn on the cooler:** Ensure that the cooler inner and outer lids are secured on the cooler. To turn the cooler on and begin cooling, place the cold key on the circular valve at the top rim of the cooler. The valve will open with an audible 'click' and the ON indicator will show on the display. The vaccine storage compartment temperature, shown on the cooler lid, should begin to decrease within 10 minutes of turning on the cooler. Allow 20-40 minutes for cool-down to complete and the vaccine storage compartment to reach a temperature between 2 and 8 °C. Leave the cold key on the valve throughout use.

NOTICE: Do not remove the cold key prior to removing vaccines. Doing so may result in damage to the vaccines.


- 3) Load vaccines:** Once the vaccine storage compartment temperature has stabilized between 2 and 8 °C, open the cooler inner lid to place vaccines in the vaccine storage compartment. Place vaccines below the "FILL LINE", indicated by the label. Close the cooler inner lid. Vaccine vials may be placed directly into the vaccine storage compartment. Secondary packaging or padding may also be used to stabilize vials during transit if space permits.

NOTICE: Do not remove outer lid while vaccines are in the vaccine storage compartment. Doing so will result in inaccurate lid sensor readings and may result in damage to vaccine vials.

NOTICE: Do not fill vaccines above the fill line. Doing so may result in damage to vaccines.

- 4) During use:** Once cool-down is complete, the vaccine carrier will maintain a temperature between 2 and 8 °C for five days of continuous use at 43 °C ambient temperature. See “Intermittent Use” section for details on non-continuous use. The vaccine storage compartment temperature and estimated remaining cold life will be shown on the cooler lid display. The cooler inner lid may be opened to access vaccines. Keep the cooler inner lid closed whenever possible during use to maintain the vaccine storage compartment temperature. The vaccine carrier may be operated and transported in any orientation. The cooler should always be transported and stored inside of the bag. Minimize shaking of the vaccine carrier to avoid damage to vaccine vials.


NOTICE: Do not keep vaccine vials in the vaccine carrier after the cold life indicator shows “0 (zero) DAYS”, as the vaccine storage compartment may not stay between 2 and 8 °C, even while the cold key is still on the valve. Leaving vaccines in the vaccine carrier may result in damage to the vaccines.

- 5) Remove any remaining vaccines:** Prior to the display showing “0 (zero) DAYS” of cold life, remove any unused vaccine vials from the vaccine storage compartment and store them in another refrigeration device or an Indigo Vaccine Carrier with remaining cold life.

- 6) After use:** After any remaining vaccines have been removed from the device, remove the cold key from the valve to turn off the cooler. Collect data from

the cooler lid using a compatible data device as required (see "Data Download"). Charge the cooler prior to next use (see "Charging").

NOTICE: The vaccine storage compartment temperature may stay between 2 and 8 °C for longer than five days at ambient temperatures lower than 43 °C.

Intermittent Use


The Indigo Vaccine Carrier is designed for five days of continuous use at 43 °C ambient temperature. The cooler may be turned off at the discretion of the user when vaccines are not present, although additional cool-downs from ambient temperature after the initial cool-down will deplete the cooler's cold life. To account for this, the cooler lid will subtract a maximum of one day from the cold life indicator each

time the vaccine carrier is turned off long enough for the vaccine chamber temperature to rise above 10 °C. Users should ensure any remaining vaccines are returned to another refrigeration device or an Indigo Vaccine Carrier with remaining cold life prior to the display showing "0 (zero) DAYS" of cold life. Upon recharging the cooler, the cold life indicator will return to its normal countdown.

Bag Strap Adjustment

To adjust the bag's straps, put the bag on your back and pull down on the free ends of the straps, near the top and bottom of the bag.

To loosen the bag straps, pull upward on the black plastic strap adjusters near the top and bottom of the straps.

Data Download

The cooler lid will store the past two years of data, which can be transferred to a compatible data device. Although any USB mass storage device can be used to download data, it is recommended to use a compatible equipment monitoring device, such as the Varo™ app. As the cooler lid stores data using relative time, absolute time will only be automatically recorded with the data when a compatible equipment monitoring device is used. The following data will be available for download from the cooler lid:

- **60-Day User-Readable Summary:** A report showing device summary information for the previous 60 days, in .pdf format.
- **60-Day Current Data:** A report showing 15-minute resolution vaccine storage compartment temperature, ambient temperature, main valve sensor, and cooler lid sensor data, in .json format.
- **Data History Folder:** Folder containing raw data for all monitored parameters for the previous two years, in .json format.

To access data using Android device with Varo™ Cold Chain Reporting App:


Use your mobile device to access Indigo data and send an email report in 2 minutes or less.

- 1) Download Varo Cold Chain Reporting App** from the Google Play store.
- 2) Open Varo:** Follow the instructions in the application. No account is needed.

- 3) Connect your phone/tablet to the cooler lid:** Remove the USB-A cable from the cooler lid USB slot and connect lid to the Android device using a USB On-The-Go (OTG) adaptor. An OTG adaptor will be included in the front pocket of the bag. If you do not have an OTG adaptor or are experiencing issues with connectivity, Varo will walk you through how to manually collect a useful snapshot report.

1


Find more information on Varo at: Varo-App.org

Battery Replacement

The battery in the cooler lid is intended to provide power for approximately five years. A spare CR123A battery is stored in the cooler lid electronics enclosure.

When the replace battery icon displays on the screen, follow these instructions to replace the battery. If the display does not turn on when the cooler is turned on with the cold key, this may also indicate a battery replacement is needed.

- 1) Using a Phillips M3 screwdriver, loosen and remove the five (5) M3 screws on the electronics enclosure.
- 2) Lift enclosure lid carefully and flip enclosure lid over to locate the connected battery.

NOTICE: Do not pull wires between the two enclosure parts. Doing so may result in damage to the lid electronics.

- 3) If this is the first time the battery is replaced, there will be a spare battery inside the enclosure. Otherwise, a new CR123A battery will need to be acquired.
- 4) Remove and safely dispose of the discharged CR123A battery and replace with a new battery. Insert batteries in the direction indicated by the plus (+) and minus (-) signs on the enclosure.
- 5) Carefully restore the enclosure lid to its original location, ensuring no wires are caught under the battery.

- 6) Return and tighten the five (5) M3 screws to the electronics enclosure lid.
- 7) Verify that the low battery icon is no longer shown on the display.

Cleaning

To clean the vaccine carrier, remove the cooler and cooler lid from the bag and empty any contents of the bag's pockets. To clean the bag, rinse with clean water. Unzip all zippers and leave in a warm, clean

environment until dry. To clean the cooler lid and cooler, use a damp towel and mild soap to wipe down all surfaces. Dry all components fully before placing back in the bag.

NOTICE: Do not remove any linings or internal foam components. Doing so may result in damage to the bag.

NOTICE: Do not use any highly acidic, highly alkaline, or flammable cleaning agents. Do not pressure-wash or use an abrasive brush to clean. Doing so may result in damage to the device.

Signs of Damage

The following may indicate a damaged cooler lid:

- Damage to the cooler outer lid that impedes its ability to be secured on the cooler.
- Damage to the cooler inner lid that impedes its ability to fit in the cooler outer lid.
- Visible damage to the electronics enclosure or display that impedes proper use, such as poor display contrast or screen cracks.
- Temperature probe has exposed wires or is loose and able to rotate.
- Nothing is shown on the cooler lid display when the cooler inner lid is removed, or the cold key is on the valve. See "Troubleshooting" section for troubleshooting steps.
- The ON indicator does not appear on the display when the cooler outer lid is on the cooler and the cold key is on the valve.

- The lid OPEN indicator does not appear on the display when the cooler inner lid is removed.
- The display shows a temperature outside the realistic range (-10 to 50 °C) or reads “Err”.
- When the cooler lid USB is connected to an Android device with the Varo App, the Varo App does not show that the cooler lid has been connected. See “Troubleshooting” section for troubleshooting steps.

The following may indicate a damaged bag:

- Any damaged components that prevent proper use, such as jammed zippers or torn straps.

The following may indicate a damaged cooler:

- Vaccine storage compartment temperature rises above 8 °C while the cooler lid display indicates remaining cold life.
- The vaccine carrier is not cooling after the cold key has been placed on the valve. See “Troubleshooting” section for troubleshooting steps.
- Cool-down takes longer than two hours.
- The cooler does not fit into the charger properly.
- Significant damage (denting, punctures, etc.) is visible on the cooler.

If any of the above occur and are not resolved by troubleshooting steps, contact

the service provider and submit a warranty claim, if applicable (see “Warranty Terms and Conditions”). See the service contact card in

the card holder on the top inside panel of the bag or visit indigocold.org/support for service contact information.

Troubleshooting

If you experience issues with the Indigo Vaccine Carrier, try the following troubleshooting steps. If the following steps are not successful, contact the service provider and submit a warranty claim, if applicable (see “Warranty Terms and Conditions”). See the service contact card in the card holder on the top inside panel of the bag or visit indigocold.org/support for service contact information.

The vaccine carrier is not cooling after the cold key has been placed on the valve.

- Adjust the cold key on the valve to make sure it is seated correctly. The valve should make an audible “click” upon placement and should not be tilted.
- Charge the cooler and retry cool-down.

Nothing is shown on the cooler lid display when the cooler inner lid is removed, or the cold key is on the valve.

- The cooler will continue work without a functioning lid, but lid features will not be available, such as temperature monitoring and cold life estimates.
- The display should turn on if the cooler inner lid is removed or the cold key is on


the valve. Check that the cooler outer lid has been firmly pressed on the cooler and that the cooler inner lid is fully removed, or the cold key is on the valve.

- Replace the cooler lid battery. See “Battery Replacement” section for instructions.

When the cooler lid USB is connected to an Android device with the Varo App, the Varo App does not show that the cooler lid has been connected.

- Check that the cooler lid is connected to the Android device using a USB On-The-Go (OTG) adaptor. A USB OTG adaptor will be included in the front pocket of the bag.
- Check cooler lid USB connection using an alternate device, either an Indigo Charger or a different Android device with Varo App. If checking with a different Android device, be sure to use a different USB OTG adaptor as well. If the lid does not register as connected with the alternate device, there may be an issue with the cooler lid. If the lid registers as connected with the alternate device, use a different USB OTG adaptor and Android device with Varo App for data download.

CHARGING

Safety Precautions

WARNING: Read all instructions carefully before use. These precautions will reduce risk of burns, electric shock, fire, and injury to persons. Follow basic safety precautions, including those in the following pages. This appliance is not intended for use by persons, including children, with reduced physical, sensory, or mental capabilities, or lack of experience and knowledge.

Proper Installation and Maintenance

- Remove all packaging prior to use.
- The charger is intended for indoor use only. It is not approved for outdoor use.
- Do not store or use corrosive chemicals, vapors, or flammables in or near the charger.
- Do not obstruct air vents.
- Do not operate the charger if it is damaged or not working properly.
- Do not repair or replace parts unless instructed to do so in this manual.


Burn Prevention

Do not open the charger until the orange “WARNING: HOT, DO NOT TOUCH” light has turned off and the green “COMPLETE” light has turned on.

- During charging, interior surfaces of the charger are hot enough to cause burns. Do not open the charger until it has had sufficient time to cool. The charger will indicate that cooling is complete by the “WARNING: HOT” light turning off and the “COMPLETE” light turning on. Cooling should take one to two hours.
- During charging, do not touch or let any flammable materials touch the interior surfaces of the charger.
- Do not leave the charger unattended during charging.
- Secure all loose garments, long hair, etc. prior to use.

Fire Safety

- Use the charger only for its intended use as described in this manual.
- Remove anything in or attached to the cooler, other than the tether and labels on the cooler rim, prior to charging.
- A working smoke detector should be near the charger while in use.
- Have an appropriate fire extinguisher available, nearby, highly visible, and easily accessible.
- Do not leave the charger unattended during the charging cycle.

Child Safety

- Do not allow children to use the charger. The charger is intended for use by trained operators and is not intended for use by children.
- Do not allow climbing, standing, leaning, sitting, or hanging on any part of the charger.
- Do not allow children in the vicinity of the charger while it is in use. Children should not be left alone or unattended

in the area when the charger is not in use. Children should never be allowed to play in the vicinity of the charger.


Cleaning Safety

- Follow cleaning and maintenance instructions in this manual.
 - Unplug device prior to cleaning.
-

Components


Overview

- 1) Charger Lid Conduit
- 2) Charger Lid Handles
- 3) Handles
- 4) Housing
- 5) Vents
- 6) (4) Legs
- 7) Charger Lid Fan
- 8) Charger Lid Holder
- 9) Cooler Lid Holder
- 10) Display
- 11) Charging and USB Ports
- 12) Cooler Chamber


Display

- A) CHARGING Light
- B) START Button
- C) COOLING Light
- D) COMPLETE Light
- E) Screen
- F) ERROR Light
- G) WARNING: HOT Light

Charging and USB Ports


- H) LID Port
- I) DATA Port
- J) (2) CHARGER 5V Ports

Installation

Shipped Components and Tools

- 1) Charger Body
- 2) (4) Legs
- 3) Phillips Screwdriver
- 4) Charger Lid Holder
- 5) Cooler Lid Holder


Installation Instructions


NOTICE: If charger is repackaged after initial installation, the lid latches should be reattached to the charger lid to prevent damage to the device during shipping.

Part 1: Preparation and Lid Latch Removal


- 1) Read instruction manual carefully.
- 2) Remove the top packaging pieces, charger legs, charger lid holder, and cooler lid holder from the shipping box.
- 3) With two people lifting the black charger lid handles, remove the charger body from the box.
- 4) Remove screw attaching the lid latch to the charger lid using the Phillips screwdriver.

- 5) Turn lid latch towards threaded hole on charger body.


- 6) Install previously removed screw through open hole on lid latch into threaded hole on charger body using Phillips screwdriver.


- 7) Repeat steps 4-6 on the lid latch on the other side of the charger lid.
- 8) Check that all components and tools are present (see "Shipped Components and Tools").

Part 2: Leg Installation


- 1) Fit the larger holes on the leg over the screw heads on the charger body. The charger may be placed on its side to install legs. Use the Phillips screwdriver to remove the screws if necessary.

- 2) Press the legs towards the top of the charger so that the screw heads sit in the smaller holes. Use the Phillips screwdriver to loosen the screws if necessary.


- 3) Tighten four (4) screws using Philips screwdriver to secure legs.


- 4) Repeat steps 1-4 for all four charger legs.

Part 3: Charger Lid Holder Installation


- 1) Fit larger holes on charger lid holder over screw heads on charger body.

- 2) Press charger lid holder towards the bottom of the charger so that the screw heads sit in the smaller holes. Use the Phillips screwdriver to loosen the screws if necessary.


- 3) Tighten four (4) screws using Phillips screwdriver to secure charger lid holder.


Part 4: Cooler Lid Holder Installation


1) Fit the larger holes on the cooler lid holder over the screw heads on the charger body.

2) Press the cooler lid holder towards the bottom of the charger so that the screw heads sit in the smaller holes. Use the Phillips screwdriver to loosen the screws if necessary.


3) Tighten two (2) screws using Phillips screwdriver to secure cooler lid holder.

Part 5: Preparation for Use

WARNING: Do not obstruct air vents. Doing so may result in fire and injury to persons.

- 1) Place the charger in a dry, covered location with at least 15 cm (6 inches) of space around vents and enough space to open the charger lid.
- 2) Open the charger lid to remove foam and cardboard packaging from charger body.
- 3) Plug power cable into a 220-240V AC power source and check that circuit breaker is in the "on" position. See "Power Requirements" section for more details.

WARNING: Do not use charger without removing foam and cardboard packaging from charger body. Failure to remove packaging may result in fire and injury to persons.

Power Requirements

NOTICE: This appliance may deplete or drain power from facility systems that were not intended or designed for this appliance.

Nominal Supply Voltage: 220 - 240 V AC

Nominal Current Draw: 5.4 - 3.9 A

Nominal Supply Frequency: 50 - 60 Hz

Phases: 1

Operational Supply Voltage: 110 - 278 V AC

Minimum Current Draw: 3.3 A (278 V AC)

Maximum Current Draw: 8.2 A (110 V AC)

Operation



NOTICE: Cooler tether can withstand charging temperatures and should remain attached to the cooler during charging.

NOTICE: Upon startup, the charger may give off an odor. This is considered normal, although signs of fire such as smoke are abnormal and should be addressed immediately.


- 1) Confirm charger is in a dry, covered location and plugged in, and circuit breaker is in the “on” position. Open charger lid and place in charger lid holder.

To open charger lid:


- A) Lift the charger lid slightly.
- B) Rotate counterclockwise until the arrows on the charger lid and charger body align.
- C) Lift the charger lid out of the charger body.

- 2) Place the vaccine carrier, including cooler, cooler lid, and bag, to the right of the charger (from the perspective of someone standing in front of the charger).


- 3) Unzip the bag's lower zipper and open the bag.


NOTICE: Do not remove cooler tether from cooler.

- 4) Remove the full cooler lid (both inner and outer lids) from cooler and place on cooler lid holder.

WARNING: Do not leave anything, other than the tether and labels on the cooler rim, in or attached to the cooler during charging. Doing so may result in fire and injury to persons.

- 5) Loosen the cinch cord from around the cooler. If necessary, use the orange tabs to pull the liner away from the cooler. Check that anything in or attached to the cooler, other than the tether and labels on the cooler rim, have been removed prior to charging.


NOTICE: Do not lift cooler by the valve. Doing so may result in damage to the cooler.

- 6) Lift the cooler by the upper rim, place in the charger, and line up valve with arrow at front of the charger.

INSERT COOLER
AND CLOSE CHARGER
LID TO BEGIN

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
NOTICE: Ensure that the cooler's tether is not caught on the charger lid tabs. Failure to do so may result in damage to the tether.

- 7) Remove the charger lid from the charger lid holder and line up the arrow on the charger lid with arrow on the charger body. Twist the charger lid clockwise to move into the charging position. If done correctly, the charger lid will drop into place so that it is flush with the charger body. The charger will not start if the charger lid is not in the charging position.

INSERT COOLER
AND CLOSE CHARGER
LID TO BEGIN

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
- 8) Remove USB-A cable from the cooler lid USB slot and plug into the charger USB-A port.

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
- 9) Press START button to begin charge. Charging a fully discharged cooler will take approximately four to five hours.

READY TO CHARGE
PRESS BUTTON
TO BEGIN

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
WARNING: Do not open charger lid or remove cooler while “WARNING: HOT” light is on. Doing so may result in injury.

- 10) When the screen shows “CHARGING COMPLETE SAFE TO TOUCH” and the “WARNING: HOT” light has turned off, it is safe to open the charger lid. To open the charger lid, lift the charger lid slightly, rotate counterclockwise until the arrows on the charger lid and charger body align, and lift the charger lid out of the charger body. Place the charger lid in the charger lid holder.

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

NOTICE: Do not place the cooler in the bag or put the cooler lid on the cooler while the cooler is hot. Doing so may result in damage to the bag or cooler lid.


- 11) Check that the bag is placed to the right of the charger, the lower zipper is unzipped, and the bag lid is open. Lift the cooler by the cooler rim and place the cooler in the bag. Tighten cinch cord around the cooler.

- 12) Unplug the cooler lid USB cable from the charger, place USB cable back into the cooler lid USB slot.


- 13) Place full cooler lid (both inner and outer lids) back on the cooler. Ensure that the cooler outer lid valve slot is aligned with the cooler valve and that the cold key is in its holder on the cooler lid (not on valve).

- 14) Open the cooler inner lid to display the cold life indicator and check the cooler's charge. A fully charged cooler will show "5 DAYS" of cold life with a complete circle around the number.

- 15) Close bag and zip shut.

- 16) Store vaccine carrier in a dry location without the cold key on the valve until use.

NOTICE: If a charged vaccine carrier is not used within six months, repeat the steps above to recharge the vaccine carrier.

Data Download

The charger will store the past two years of data, which can be transferred to a compatible data device. Although any USB mass storage device can be used to download data, it is recommended to use a compatible equipment monitoring device, such as the Varo™ app. As the charger stores data using relative time, absolute time will only be automatically recorded with the data when a compatible

equipment monitoring device is used. The following data will be available for download from the charger:

- **60-Day Current Data:** A report showing raw data for all monitored parameters for the previous 60 days, in .json format.
- **Data History Folder:** Folder containing raw data for all monitored parameters for the previous two years, in .json format.

To access data using Android device with Varo™ Cold Chain Reporting App:

Use your mobile device to access Indigo data and send an email report in 2 minutes or less.


- 1) Download Varo Cold Chain Reporting App** from the Google Play store.
- 2) Open Varo:** Follow the instructions in the application. No account is needed.
- 3) Connect your phone/tablet to the charger:** Connect the Android device to the charger using a USB On-The-Go

(OTG) adaptor, USB-A to Micro-USB cable (USB cable), and the charger Micro-USB "DATA" port. Plug the Micro-USB side of the USB cable into the charger Micro-USB "DATA" port. Connect the USB-A side of the USB cable to the Android device using the OTG adaptor. An OTG adaptor will be included with the charger. If you do not have a OTG adaptor, USB cable, or are experiencing issues with connectivity, Varo will walk you through how to manually collect a useful snapshot report.

Scan QR Code
Varo-App.org

1

2

3

Find more information on Varo at:
Varo-App.org

Cleaning

If the charger is dirty, disconnect from power, then use a damp towel to wipe down exposed surfaces including the charger lid fan cover and electronics enclosure vents.

NOTICE: Do not use any highly acidic, highly alkaline, or flammable cleaning agents. Do not pressure-wash or use an abrasive brush to clean. Doing so may result in damage to the device.

Signs of Damage

The following may indicate a damaged Charger:

- Charger screen displays “ERROR” or “FAULT” and/or the red “ERROR” light is on.
- Nothing is shown on the charger display while it is connected to power.
- Visual damage to the charger display that impedes proper use, such as poor display contrast or screen cracks.
- Damage to the power cord.

NOTICE: If the supply cord is damaged, it must be replaced by the manufacturer, its service agent, or similarly qualified persons in order to avoid a hazard.

- Damage to the charger such that the cooler does not fit into the charger.
- Damage causing exposed insulation or loss of insulation.
- Damage to the charger lid such that it does not fit on the charger or charger does not move past “CLOSE CHARGER LID TO BEGIN” screen when charger lid is placed on the charger.
- See “Troubleshooting” section for troubleshooting steps.
- When the charger USB is connected to an Android device with the Varo App, the Varo App does not show that the charger has been connected.
- When the cooler lid USB is connected to the charger, the charger does not move past the “CONNECT COOLER LID USB

LID USB TO BEGIN” screen. See “Troubleshooting” section for troubleshooting steps.

- Cooler lid cold life indicator does not update to “5 DAYS” after cooler lid is connected during charging and charger display shows “CHARGE COMPLETE”.

If any of the above occur and are not resolved by troubleshooting steps, contact the service provider and submit a warranty claim, if applicable (see “Warranty Terms and Conditions”). See the service contact card in the card holder on the top inside panel of the bag or visit indigocold.org/support for service contact information.

Troubleshooting

If you experience issues with the Indigo Charger, try the following troubleshooting steps. If the following steps are not successful, contact the service provider and submit a warranty claim, if applicable (see “Warranty Terms and Conditions”). See the service contact card in the card holder on the top inside panel of the bag or visit indigocold.org/support for service contact information.

Charger does not move past “CLOSE CHARGER LID TO BEGIN” screen when charger lid is placed on the charger.

- Check that the charger lid is fully seated on the charger. If done correctly, the charger lid will drop into place so that it is flush with the charger body. The charger will not start if the charger lid is not in the charging position.

When the cooler lid USB is connected to the charger, the charger does not move past the “CONNECT COOLER LID USB TO BEGIN” screen.

- First, determine if the issue lies with the USB connection on the cooler lid. Connect the cooler lid to an Android device with the Varo app using a USB OTG adaptor. An OTG adaptor will be

included in the front pocket of the bag. If the lid does not register as connected with the device, the issue is with the cooler lid USB connection.

- If the lid registers as connected, next, determine if the issue lies with the USB connection on the charger. Connect a different cooler lid to the charger. If the cooler lid does not register as connected with the charger, the issue is with the charger USB connection.
- In either of these USB issue scenarios, you may continue to use the vaccine carrier and/or charger for immunization activities, but lid cold life will not update, and data collection will not be possible. Please note that charging will not begin automatically when cooler is loaded to the charger in this scenario. To manually initiate charging, load the cooler into the charger per the usual procedure. When the charger screen displays “CONNECT COOLER LID USB TO BEGIN”, press and hold the “START” button on the charger for five seconds. This will initiate the charging process.
- If both devices can connect to an Android device with the Varo app, the issue is likely with the charger.

DISPOSAL AND RECYCLING GUIDELINES

When the Indigo Vaccine Carrier is no longer in use, to fully decommission, remove the cooler from the bag and cut the tether from the cooler lid to the cooler. Label the cooler as “Non-functioning” before discarding.

When an Indigo Charger is no longer in use, decommission it by unplugging it and labeling it “Non-functioning”.

WARNING: Do not attempt to remove or modify the insulation within the charger heating enclosure. Doing so may result in irritation of the skin, eyes, and lungs.

WARNING: Do not attempt to remove or modify the electronics within the charger, as this is a high voltage electrical hazard. Doing so may result in electric shock.

Disposal of the Indigo system must be made according to local appliance disposal and recycling regulations.

Major materials in the vaccine carrier include:

- **Stainless steel (cooler)**
- **Desiccant (cooler)**
- **Water (cooler)**
- **Lithium-ion battery (cooler lid)**

- **Polypropylene foam (cooler lid)**
- **Polycarbonate (cooler lid)**
- **Nylon fabric (bag)**
- **Polyethylene foam (bag)**

Major materials in the charger include:

- **Stainless steel**
- **Fiberglass insulation**
- **Mica insulation**

Contact your service provider for further information. See the service contact card in the card holder on the top inside panel of the bag or visit indigocold.org/support for service contact information.

WARRANTY TERMS AND CONDITIONS

Contact your service provider to make a warranty claim. See the service contact card in the card holder on the top inside panel of the bag or visit indigocold.org/support for service contact information. The following coverage will be provided:

Two-year warranty on the Indigo Vaccine Carrier and Charger, beginning at the date of shipment from manufacturer, included at no extra cost. Under warranty coverage, the manufacturer will repair or replace device components that have failed in the use of the device for its designed intent.

Free warranty coverage does not include the following cases:

- Damage caused by improper use or improper maintenance.
 - Damage caused by attempted repairs done by a non-designated repair department.
 - Damage caused by force majeure.
 - Damage occurring after the warranty period has expired.
-

ICON GLOSSARY

Icon	Icon Descriptions
	Cold key activates valve through magnetic force.
	Do not fill vaccines above the fill line.
	Do not pull wires between the two enclosure parts.
	Hot, do not touch.

Icon	Icon Descriptions
	Do not obstruct air vents.
	Do not lift the cooler by the valve.
	Align valve with arrow at front of the charger.
	USB-A cable

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